

RAB2B Antibody



PACO05546

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

RAB2B

Uniprot

Q8WUD1

Synonyms:

RAB2B; Ras-related protein Rab-2B

Immunogen:

Synthesized peptide derived from the Internal region of human Rab 2B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A