

PACO05552

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

PIGY

Uniprot

Q3MUY2

Synonyms:

PIGY; Phosphatidylinositol N-acetylglucosaminyltransferase subunit Y;
Phosphatidylinositol-glycan biosynthesis class Y protein; PIG-Y

Immunogen:

Synthesized peptide derived from the N-terminal region of human PIG-Y.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A