

REPS1 Antibody



PACO05557

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

REPS1

Uniprot

Q96D71

Synonyms:

REPS1; RalBP1-associated Eps domain-containing protein 1; RalBP1-interacting protein 1

Immunogen:

Synthesized peptide derived from the C-terminal region of human REPS1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A