

# BAGE4 Antibody



PACO05564

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

BAGE4

**Uniprot**

Q86Y28

**Synonyms:**

BAGE4; MLL3P; B melanoma antigen 4; Cancer/testis antigen 2.4; CT2.4

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human BAGE4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A