

Phospho-CAV2 (Y27) Antibody



PACO05582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

CAV2

Uniprot

P51636

Synonyms:

CAV2; Caveolin-2

Immunogen:

Synthesized peptide derived from human Caveolin-2 around the phosphorylation site of Y27.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A