

SERPINA6 Antibody



PACO05596

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

SERPINA6

Uniprot

P08185

Synonyms:

SERPINA6; CBG; Corticosteroid-binding globulin; CBG; Serpin A6; Transcortin

Immunogen:

Synthesized peptide derived from the Internal region of human SerpinA6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A