

# PEA15 Antibody



PACO05606

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PEA15

**Uniprot**

Q15121

**Synonyms:**

PEA15; Astrocytic phosphoprotein PEA-15; 15 kDa phosphoprotein enriched in astrocytes; Phosphoprotein enriched in diabetes; PED

**Immunogen:**

Synthesized peptide derived from human PEA-15 around the non-phosphorylation site of S104.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A