

B4GALT3 Antibody



PACO05609

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

B4GALT3

Uniprot

O60512

Synonyms:

B4GALT3; beta-1; 4-galactosyltransferase 3; beta-1,4-GalTase 3; beta4Gal-T3; b4Gal-T3;
UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 3; UDP-galactose: beta-N-
acetylglucosamine beta-1,4-galactosyltransferase 3

Immunogen:

Synthesized peptide derived from the Internal region of human β-1,4-Gal-T3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A