

RIOK3 Antibody



PACO05621

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

RIOK3

Uniprot

O14730

Synonyms:

RIOK3; SUDD; Serine/threonine-protein kinase RIO3; RIO kinase 3; sudD homolog

Immunogen:

Synthesized peptide derived from the Internal region of human SUDD.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A