

Phospho-CCND2 (T280) Antibody



PACO05648

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

CCND2

Uniprot

P30279

Synonyms:

CCND2; G1/S-specific cyclin-D2

Immunogen:

Synthesized peptide derived from human Cyclin D2 around the phosphorylation site of T280.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A