

Phospho-CCND3 (T283) Antibody



PACO05650

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

CCND3

Uniprot

P30281

Synonyms:

CCND3; G1/S-specific cyclin-D3

Immunogen:

Synthesized peptide derived from human Cyclin D3 around the phosphorylation site of T283.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A