

# Phospho-CCNH (T315) Antibody



PACO05665

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

CCNH

**Uniprot**

P51946

**Synonyms:**

CCNH; Cyclin-H; MO15-associated protein; p34; p37

**Immunogen:**

Synthesized peptide derived from human Cyclin H around the phosphorylation site of T315.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A