

TAAR5 Antibody



PACO05669

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

TAAR5

Uniprot

O14804

Synonyms:

TAAR5; PNR; Trace amine-associated receptor 5; TaR-5; Trace amine receptor 5; Putative neurotransmitter receptor

Immunogen:

Synthesized peptide derived from the C-terminal region of human TAAR5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A