

# DOK2 Antibody



PACO05672

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

DOK2

**Uniprot**

O60496

**Synonyms:**

DOK2; Docking protein 2; Downstream of tyrosine kinase 2; p56(dok-2)

**Immunogen:**

Synthesized peptide derived from human Dok-2 around the non-phosphorylation site of Y345.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A