

MAP3K6 Antibody



PACO05678

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

MAP3K6

Uniprot

O95382

Synonyms:

MAP3K6; ASK2; MAPKKK6; MEKK6; Mitogen-activated protein kinase kinase kinase 6; Apoptosis signal-regulating kinase 2

Immunogen:

Synthesized peptide derived from the Internal region of human MEK Kinase-6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A