

# Phospho-ZNF598 (Y306) Antibody



PACO05684

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

ZNF598

**Uniprot**

Q86UK7

**Synonyms:**

ZNF598; Zinc finger protein 598

**Immunogen:**

Synthesized peptide derived from human ZFP598 around the phosphorylation site of Y306.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A