

TBX18 Antibody



PACO05688

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

TBX18

Uniprot

O95935

Synonyms:

TBX18; T-box transcription factor TBX18; T-box protein 18

Immunogen:

Synthesized peptide derived from the Internal region of human TBX18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A