

# Phospho-LATS1/LATS2 (T1079/1041) Antibody



PACO05691

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

LATS1/LATS2

**Uniprot**

O95835/Q9NRM7

**Synonyms:**

LATS1; WARTS; Serine/threonine-protein kinase LATS1; Large tumor suppressor homolog 1; WARTS protein kinase; h-warts; LATS2; KPM; Serine/threonine-protein kinase LATS2; Kinase phosphorylated during mitosis protein; Large tumor suppressor ho

**Immunogen:**

Synthesized peptide derived from human LATS1/2 around the phosphorylation site of T1079/1041.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A