

# BTF3L4 Antibody



PACO05697

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

BTF3L4

**Uniprot**

Q96K17

**Synonyms:**

BTF3L4; Transcription factor BTF3 homolog 4; Basic transcription factor 3-like 4

**Immunogen:**

Synthesized peptide derived from the Internal region of human BTF3L4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A