

COL23A1 Antibody



PACO05701

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:10000, IF:1:200-1:1000

Protein Background:**Gene ID:**

COL23A1

Uniprot

Q86Y22

Synonyms:

COL23A1; Collagen alpha-1(XXIII) chain

Immunogen:

Synthesized peptide derived from the C-terminal region of human COL23A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A