

Phospho-CD4 (S433) Antibody



PACO05717

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CD4

Uniprot

P01730

Synonyms:

CD4; T-cell surface glycoprotein CD4; T-cell surface antigen T4/Leu-3; CD antigen CD4

Immunogen:

Synthesized peptide derived from human CD4 around the phosphorylation site of S433.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A