

CD6 Antibody



PACO05724

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CD6

Uniprot

P30203

Synonyms:

CD6; T-cell differentiation antigen CD6; T12; TP120; CD antigen CD6

Immunogen:

Synthesized peptide derived from the Internal region of human CD6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A