

# PEX11G Antibody



PACO05740

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

PEX11G

**Uniprot**

Q96HA9

**Synonyms:**

PEX11G; PEX11C; Peroxisomal membrane protein 11C; Peroxin-11C; Peroxisomal biogenesis factor 11C; Protein PEX11 homolog gamma; PEX11-gamma

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human Peroxin 11 &gamma; .

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A