

GTF3C4 Antibody



PACO05746

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

GTF3C4

Uniprot

Q9UKN8

Synonyms:

GTF3C4; General transcription factor 3C polypeptide 4; TF3C-delta; Transcription factor IIIC 90 kDa subunit; TFIIC 90 kDa subunit; TFIIC90; Transcription factor IIIC subunit delta

Immunogen:

Synthesized peptide derived from the C-terminal region of human TFIIC90.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A