

NRXN3 Antibody



PACO05751

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

NRXN3

Uniprot

Q9HDB5

Synonyms:

NRXN3; KIAA0743; Neurexin-3-beta; Neurexin III-beta

Immunogen:

Synthesized peptide derived from the N-terminal region of human Neurexin III β .

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A