

NRXN1 Antibody



PACO05753

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

NRXN1

Uniprot

Q9ULB1

Synonyms:

NRXN1; KIAA0578; Neurexin-1-alpha; Neurexin I-alpha

Immunogen:

Synthesized peptide derived from the Internal region of human Neurexin I.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A