

# ARHGAP12 Antibody



PACO05758

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ARHGAP12

**Uniprot**

Q8IWW6

**Synonyms:**

ARHGAP12; Rho GTPase-activating protein 12; Rho-type GTPase-activating protein 12

**Immunogen:**

Synthesized peptide derived from the Internal region of human ARHGAP12.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A