

MRPL33 Antibody



PACO05777

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:

Gene ID:

MRPL33

Uniprot

O75394

Synonyms:

MRPL33; C2orf1; 39S ribosomal protein L33; mitochondrial; L33mt; MRP-L33

Immunogen:

Synthesized peptide derived from the N-terminal region of human MRP-L33.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A