

# ZNF592 Antibody



PACO05796

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

ZNF592

**Uniprot**

Q92610

**Synonyms:**

ZNF592; KIAA0211; Zinc finger protein 592

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human ZNF592.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A