

Phospho-MDC1 (S513) Antibody



PACO05798

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

MDC1

Uniprot

Q14676

Synonyms:

MDC1; KIAA0170; NFBD1; Mediator of DNA damage checkpoint protein 1; Nuclear factor with BRCT domains 1

Immunogen:

Synthesized peptide derived from human MDC1 around the phosphorylation site of S513.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A