

# RAB11FIP3 Antibody



PACO05806

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

RAB11FIP3

**Uniprot**

O75154

**Synonyms:**

RAB11FIP3; ARFO1; KIAA0665; Rab11 family-interacting protein 3; FIP3-Rab11; Rab11-FIP3; Arfophilin-1; EF hands-containing Rab-interacting protein; Eferin; MU-MB-17.148

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Rab11-FIP3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A