

Phospho-ACAP1 (S554) Antibody



PACO05811

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ACAP1

Uniprot

Q15027

Synonyms:

ACAP1; CENTB1; KIAA0050; Arf-GAP with coiled-coil; ANK repeat and PH domain-containing protein 1; Centaurin-beta-1; Cnt-b1

Immunogen:

Synthesized peptide derived from human ACAP1 around the phosphorylation site of S554.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A