

Cleaved-CD97 (L530) Antibody



PACO05812

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CD97

Uniprot

P48960

Synonyms:

CD97; CD97 antigen; Leukocyte antigen CD97; CD antigen CD97

Immunogen:

Synthesized peptide derived from the Internal region of human CD97 α .

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A