

ARHGAP11A Antibody



PACO05817

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

ARHGAP11A

Uniprot

Q6P4F7

Synonyms:

ARHGAP11A; KIAA0013; Rho GTPase-activating protein 11A; Rho-type GTPase-activating protein 11A

Immunogen:

Synthesized peptide derived from the Internal region of human ARHGAP11A.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A