

Phospho-GAB2 (Y643) Antibody



PACO05824

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

GAB2

Uniprot

Q9UQC2

Synonyms:

GAB2; KIAA0571; GRB2-associated-binding protein 2; GRB2-associated binder 2;
Growth factor receptor bound protein 2-associated protein 2; pp100

Immunogen:

Synthesized peptide derived from human Gab 2 around the phosphorylation site of Y643.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A