

TLK1 Antibody



PACO05832

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:5000, IF:1:200-1:1000

Protein Background:**Gene ID:**

TLK1

Uniprot

Q9Y4F6

Synonyms:

TLK1; KIAA0137; Serine/threonine-protein kinase tousled-like 1; PKU-beta; Tousled-like kinase 1

Immunogen:

Synthesized peptide derived from human TLK1 around the non-phosphorylation site of S764.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A