

URB1 Antibody



PACO05833

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

URB1

Uniprot

O60287

Synonyms:

URB1; C21orf108; KIAA0539; NOP254; NPA1; Nucleolar pre-ribosomal-associated protein 1; Nucleolar protein 254 kDa; URB1 ribosome biogenesis 1 homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human NPA1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A