

# SEC16A Antibody



PACO05840

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

SEC16A

**Uniprot**

O15027

**Synonyms:**

SEC16A; KIAA0310; SEC16; SEC16L; Protein transport protein Sec16A; SEC16 homolog A

**Immunogen:**

Synthesized peptide derived from the Internal region of human Sec16A.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A