

CDC25A Antibody



PACO05843

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

CDC25A

Uniprot

P30304

Synonyms:

CDC25A; M-phase inducer phosphatase 1; Dual specificity phosphatase Cdc25A

Immunogen:

Synthesized peptide derived from human Cdc25A around the non-phosphorylation site of T507.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A