

PACO05856

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

CDC25C

Uniprot

P30307

Synonyms:

CDC25C; M-phase inducer phosphatase 3; Dual specificity phosphatase Cdc25C

Immunogen:

Synthesized peptide derived from human Cdc25C around the non-phosphorylation site of T48.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A