

# IGFN1 Antibody



PACO05875

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat, Pg

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:50-1:200

**Protein Background:****Gene ID:**

IGFN1

**Uniprot**

Q86VF2

**Synonyms:**

eEF1A2 binding protein

**Immunogen:**

Recombinant Protein.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---

N/A

N/A