

# ESR1 Antibody



PACO05879

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

ESR1

**Uniprot**

P03372

**Synonyms:**

Estrogen Receptor alpha

**Immunogen:**

Recombinant Protein.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---

N/A

N/A