

NFE2L2 Antibody



PACO05882

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:50-1:200

Protein Background:**Gene ID:**

NFE2L2

Uniprot

Q16236

Synonyms:

Nrf2

Immunogen:

Recombinant Protein.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A