

Bombesin Receptor 2 Antibody



PACO05886

Product Information

Size:

50ug

Reactivity:

Rat, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

Uniprot

N/A

Synonyms:

Bombesin Receptor 2

Immunogen:

Synthetic peptide.

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images

N/A

N/A