

# OR8B4 Antibody



PACO05894

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR8B4

**Uniprot**

Q96RC9

**Synonyms:**

OR8B4; OR8B4P; Olfactory receptor 8B4; Olfactory receptor OR11-315

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 8B4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A