

Phospho-IRAK1 (T100) Antibody



PACO05905

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

IRAK1

Uniprot

P51617

Synonyms:

IRAK1; IRAK; Interleukin-1 receptor-associated kinase 1; IRAK-1

Immunogen:

Synthesized peptide derived from human IRAK-1 around the phosphorylation site of T100.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A