

Phospho-IRF3 (S396) Antibody



PACO05910

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

IRF3

Uniprot

Q14653

Synonyms:

IRF3; Interferon regulatory factor 3; IRF-3

Immunogen:

Synthesized peptide derived from human IRF-3 around the phosphorylation site of S396.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A