

Phospho-AR (S650) Antibody



PACO05919

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

AR

Uniprot

P10275

Synonyms:

AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4

Immunogen:

Synthesized peptide derived from human AR around the phosphorylation site of S650.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A