

# ITPR1 Antibody



PACO05942

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

ITPR1

**Uniprot**

Q14643

**Synonyms:**

ITPR1; INSP3R1; Inositol 1; 4,5-trisphosphate receptor type 1; IP3 receptor isoform 1; IP3R 1; InsP3R1; Type 1 inositol 1,4,5-trisphosphate receptor; Type 1 InsP3 receptor

**Immunogen:**

Synthesized peptide derived from human IP3R-I around the non-phosphorylation site of S1598.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A