

# Phospho-JAK2 (Y119) Antibody



PACO05946

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

JAK2

**Uniprot**

O60674

**Synonyms:**

JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2

**Immunogen:**

Synthesized peptide derived from human JAK2 around the phosphorylation site of Y119.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A